

Ventyv[®]

Nitrile PF Exam Aloe & Vitamin E
ULTRA 3.5 Coral



Ventyv[®]
Outsmart Infection[™]



sritrangusa.com | ventyv.com | Hello@sritrangusa.com | 1.888.483.6898 |
+1 (813) 606-4301 | 5820 West Cypress Street, Suite H | Tampa, FL 33607

Ventyv[®]

Nitrile PF Exam Aloe & Vitamin E
ULTRA 3.5 Coral



Ventyv[®] Nitrile Powder-Free Ultra 3.5 Coral (Flamingo) Exam Glove is ideal for use in a large variety of clinical settings. This high end glove is a must-have for its Aloe & Vitamin E coating, which provides protection against dry and cracking skin.

Flamingo



Reorder #	Size	Gloves/ Box ¹	Boxes/ Case
10333106	XS	200	10
10334106	S	200	10
10335106	M	200	10
10336106	L	200	10
10337106	XL	180	10

PRODUCT ATTRIBUTES

- Coral
- Nitrile²
- Powder-free
- Polymer Coated
- Coated with Aloe & Vitamin E
- Fingertip Textured

THICKNESS (Single Wall)

- Finger: 0.05 mm / 1.96 mils
- Palm: 0.05 mm / 1.96 mils

PHYSICAL ATTRIBUTES³

- Length: 220 mm (XS, S)
230 mm (M, L, XL)
- Width: 70 ± 10 mm (XS)
80 ± 10 mm (S)
95 ± 10 mm (M)
110 ± 10 mm (L)
120 ± 10 mm (XL)

PHYSICAL PROPERTIES

- Tensile strength before aging: 14 MPa
- Tensile strength after aging: 14 MPa
- Elongation before aging: 500%
- Elongation after aging: 400%



80%

Inner boxes made with 80% recycled material

90%

Outer cartons made with 90% recycled material

1. By weight
2. This product is not made with Natural Rubber Latex
3. Physical Attributes meet or exceed ASTM standards for exam gloves

Ventyv[®]
Outsmart Infection[®]



sritrangusa.com | ventyv.com | Hello@sritrangusa.com | 1.888.483.6898 | +1 (813) 606-4301 | 5820 West Cypress Street, Suite H | Tampa, FL 33607