

LED RETRO-FIT UPGRADE

"Go Green" by
Upgrading your existing
Pelton LFII, Adec 6300 or
Belmont Clesta halogen light
with one of Flight Dental's
Energy Efficient LED Light-heads.



www.flightdentalsystems.com/retro

UPGRADE

- This Retro-Fit Kit installs on your existing arms and mounts and allows you to use your existing wiring and power supply.
- Easy and simple step by step instructions to upgrade your Halogen Light LED.
- Most installs can be performed in less than
 45 minute by your Local Service Technician.
- Entire Retro-Fit Kit and cost of installation is a fraction of alternate lighting solutions.

OUR FEATURES

- No need to replace eny expensive reflectors and fans
- Brighter LED Light to replace dull Halogen Lights
- Energy efficient Light Bulbs with 10,000 hours expectancy
- Control Light color with 3 Programmed Settings (Cure Safe, Surgical & Color Matching Mode)*
- Eliminate heat & bulb burn out

*IRIS LED

CALL YOUR LOCAL DEALER

21 Kenview Blvd, Brampton On L6T 5G7, CA Tel: 905.799.0517 Fax: 905.793.2406 info@flightdentalsystems.com www.flightdentalsystems.com

TORCH LED RETRO-FIT KIT UPGRADE

Part number	Description
RFL-305FA	Torch LED Retrofit Light Kits for A-dec 6300
RFL-305FPY	Torch LED Retrofit Light Kits for Pelton LFII, LFIII (Yoke Style)
RFL-305FPP	Torch LED Retrofit Light Kits for Pelton LFII, LFIII (Pin Style)
RFL-305FDCI	Torch LED Retrofit Light Kits for DCI
RFL-305FMA	Torch LED Retrofit Light Kits for MARUS
RFL-305FMI	Torch LED Retrofit Light Kits for MIDMARK



Specifications	
Color Temperature:	5500k
Light Intensity:	30,000 – 8,000 Lux
LED:	4 Direct LED's

IRIS LED RETRO-FIT KIT UPGRADE

Part number	Description
RFL-306GA	IRIS LED Retrofit Light Kits for A-dec 6300
RFL-306GPY	IRIS LED Retrofit Light Kits for Pelton LFII, LFIII (Yoke Style)
RFL-306GPP	IRIS LED Retrofit Light Kits for Pelton LFII, LFIII (Pin Style)
RFL-306GDCI	IRIS LED Retrofit Light Kits for DCI
RFL-306GMA	IRIS LED Retrofit Light Kits for MARUS
RFL-306GMI	IRIS LED Retrofit Light Kits for MIDMARK



Specifications		
Color Temperature:	4200-6000k	
Light Intensity:	35,000 – 8,000 Lux	
LED:	8 Indirect LED's	

Flight helps you UPGRADE to LED Lighting to Save Energy, Reduce Heat Discomfort and Improve your Working Light

www.flightdentalsystems.com/retro

CALL YOUR LOCAL DEALER

21 Kenview Blvd, Brampton On L6T 5G7, CA Tel: 905.799.0517 info@flightdentalsystems.com Fax: 905.793.2406 www.flightdentalsystems.com